


# OPERATING MANUAL

## Infrared Hot Plate / Stirrer

SLR



**WIGGENS**  
THE MAGIC MOTION

**WIGGENS GmbH**

Gässlesweg 22-24, 75334 Straubenhardt,  
Germany

Tel.: 0049 7248 4529088

**Wiggins Co., Ltd.**

Room 426, Hall A, Office Building M8, No.1  
Jiuxianqiao East Road, Chaoyang District,  
Beijing 100015, China

Tel: +86 400-809-2068

Fax: +86 400-809-2068-112

info@wiggins.com

www.wiggins.com

## Congratulations!

You have made an excellent choice.

WIGGENS thanks you for the trust you have placed in us.

This operating manual has been designed to help you gain an understanding of the operation and possible applications of our instruments. For optimal utilization of all functions, we recommend that you thoroughly study this manual prior to beginning operation.

## Unpacking and Inspecting

Please unpack the device carefully. Check that the package is right-side-up and then open it. Check that model of the product is one that you ordered. Check that there is no damage. If there is any damage, file a damage claim with the carrier. In the case of any damage a damage report should be requested immediately. These instructions must be followed fully for us to guarantee our full support of your claim for protecting against loss from concealed damage. The form required for filing such a claim will be provided by the carrier.

Changes without prior notification reserved

**Important:** keep operating manual for future use

# Content

1. Intended Use.....	4
2. Operator Responsibility.....	4
2.1 Disposal.....	5
2.2 CE Conformity.....	5
2.3. Technical Specifications.....	6
3. Safety Instructions.....	7
3.1. Explanation of Safety Notes.....	7
3.2. For your protection.....	8
3.3. For protection of the equipment.....	10
4. Operating Procedures.....	10
4.1. Environmental Operating Conditions.....	10
4.2. Installation.....	11
4.2.1. Installing the infrared Hot Plate/Stirrer.....	11
4.2.2. The dimensions of the infrared hot plates/stirrer (after connect PT100Teperature sensor)....	12
4.3. Operation.....	13
4.3.1. Overview of the Infrared hot plate / stirrer.....	13
4.3.2. Indicators and Functional Elements.....	14
4.3.3. Operation of the Infrared hot plate / stirrer without a Temperature Sensor.....	15
4.3.4. Operation of the Infrared hot plate / stirrer with a Temperature Sensor.....	18
4.3.5. Pt100 Temperature Sensor Calibration.....	20
4.4. The USB interface.....	22
5. Cleaning and Maintenance.....	24
5.1. Routine Cleaning.....	24
5.2. Maintenance.....	25
6. Transport and Storage.....	25
7. Accessories and Spare Parts.....	26
7.1 Temperature Sensor and Holder.....	26
7.2. Stir Bar and Retriever Examples.....	27
8. Service.....	29
8.1. Trouble-Shooting.....	29
8.2. Warranty.....	29
8.3. Contact /Technical Service.....	30

## 1. Intended Use

The SLR is a high power Infrared hot plate / stirrer used for mixing solutions and controlling their temperature in a precise, even and stable manner. The SLR features an excellent transmittance of heat energy due to the employment of an infrared heating element. It features a high power heating source for quick heat up and excellent temperature stability. The smooth and corrosion-resistant glass ceramic top plate, which has excellent thermal conductivity characteristics, can resist up to 700°C thermal shocks. The SLR also features a memory function for the stirring speed and temperature settings, which helps to enhance the effectiveness of repetitive processes. The large LCD screen displays the set and actual temperature and speed. The SLR can be connected to the included PT100 temperature sensor for precise external temperature control.

## 2. Operator Responsibility

### Use

- For mixing and/or heating liquids

### Range of use

- Laboratories
- Schools
- Pharmacies

This device is suitable for use in all areas except:

- Residential areas
- Areas that are connected directly to a low-voltage supply network that also supplies residential areas.

The safety of the user cannot be guaranteed

- if the appliance is operated with accessories that are not supplied or recommended by the manufacturer or if the appliance is operated improperly contrary to the manufacturer's specifications.

The products of *WIGGENS* ensure safe operation when installed, operated, and maintained according to common safety regulations. This section explains the potential dangers that may arise when operating the instrument and also specifies the most important safety precautions to preclude these dangers as far as possible.

- The operator is responsible for the qualification of the personnel operating the instrument.
- The personnel operating the instrument should be regularly instructed about the dangers involved with their job activities as well as measures to avert these dangers.
- Make sure all persons tasked with operating, installing, and maintaining the instrument have read and understand the safety information and operating instructions.
- When using hazardous materials or materials that could become hazardous, the instrument must be operated only by persons who are absolutely familiar with these materials and the instrument. These persons must be fully aware of possible risks.
- Only qualified personnel are authorized to perform configuration, installation, maintenance and repairs of the instrument.
- Routine operation can also be carried out by untrained personnel who should however be instructed by trained personnel.

If you have any questions concerning the operation of your instrument or the information in this manual, please contact us.

## 2.1 Disposal



At the end of its service life the instrument is to be disposed of in accordance with the local regulations specified for the disposal of electronic industry waste in an environmentally friendly manner.

## 2.2 CE Conformity



The products described in the operating instructions conform to the requirements of the following European guidelines:

Low voltage regulations with respect to legal harmonization of the member countries concerning electric devices for use within certain voltage limits.

EMC guideline with respect to legal harmonization of the member countries concerning electromagnetic compatibility.

**APPROVALS**  
European

EN61326-1: 2013, 2014/30/EU  
EN61010-1: 2010, 2014/35/EU  
EN50581: 2012, 2011/65/EU

## 2.3. Technical Specifications

Model	SLR
Mains Requirements	230V, 50Hz
Display Mode	LCD
Heating Capacity	1.8 kw
Max. Hot Plate Temperature	550°C
Min. Time to Boiling Point for 1L H <sub>2</sub> O	10 min
Speed Range	50-1500 rpm
Top Plate Area	280*280 mm
Heating Zone	Φ190 mm
Temperature Sensor Connector	PT100
Top Plate Material	Glass Ceramic
Dimensions	450* 282* 115 mm
Weight	4,1 kg
Mixing Capacity	25 L
Protection Category	IP20
Order No. for EU/UK-Plug	285416373

All measurements have been carried out at the stated voltage, frequency, and an ambient temperature of 25°C.  
 Technical changes without prior notification reserved.






WIGGENS Order Numbers consist of the Basic Order Number (BON) and the Order Number Addition (ONA) which explains different characteristics of the product that can vary from country to country. Order Numbers as stated on the product label and box label are stated as Full Order Numbers (FON), consisting of the BON followed by the ONA. For a full explanation of the ONA of your product, please ask your local WIGGENS support or refer to the Order Number Guide in the *WIGGENS* General Catalog.





## 3. Safety Instructions

### 3.1. Explanation of Safety Notes



In addition to the safety warnings listed, warnings are posted throughout the operating manual. These warnings are designated by an exclamation mark inside an equilateral triangle. "Warning of a dangerous situation (Attention! Please follow the documentation)."

Symbol	Additional term / Description
<b>Warning signs</b> 	<p>The danger is classified using a signal word. Read and follow these important instructions for averting dangers.</p> <p><b>Warning!</b>            Describes a <b>possibly</b> highly dangerous situation. If these instructions are not followed, serious injury and danger to life could result.</p> <p><b>Caution!</b>            Describes a <b>possibly</b> dangerous situation. If this is not avoided, slight or minor injuries could result. A warning of possible property damage may also be contained in the text.</p> <p><b>Notice!</b>            Describes a <b>possibly</b> harmful situation. If this is not avoided, the product or anything in its surroundings can be damaged.</p>
	<p><b>Note!</b>            Draws attention to something special.</p>
	<p><b>Important!</b>            Indicates usage tips and other useful information.</p>

## 3.2. For your protection

- Make sure you read and understand all instructions and safety precautions listed in this manual before installing or operating your instrument.
- Keep the operation instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Socket must be earthed (protective ground contact).
- Make sure the product is checked for proper condition regularly (depending on the conditions of use). Regularly check (at least every 2 months) the proper condition of the mandatory, warning, prohibition and safety labels.
- Connect the instrument to a power socket with earthing contact (PE-protective earth).
- The power supply plug serves as a safe disconnecting device from the line and must always be easily accessible.
- Do not stay in the area below the instrument.
- Never operate damaged equipment.
- Never operate instruments with damaged mains power cables.
- Observe all warning labels.
- Never remove warning labels.
- Be aware of tripping. Never route the connection cable in highly frequented areas.
- Be aware of possible cable damage. Repairs are to be carried out only by qualified service personnel
- Always turn off the instrument and disconnect the mains cable from the power source before performing any service or maintenance procedures, or before moving the instrument.
-  **Warning!** This is not an explosion proof instrument. Do not use with any highly flammable or explosive materials.
-  **Warning!** Risk of burns!
- Exercise caution when touching the housing parts and the heating plate. The heating plate can reach temperatures in excess of 500 °C. Pay attention to the residual heat after switching off. • Please make sure that the mains cable does not contact the heating plate.
-  **Warning!** Effects of the magnetic field have to be taken into account.(e.g. data storage media, cardiac pacemaker...)
- The Infrared hot plate / stirrer must only be operated in the presence of an exhaust system! (When heating with volatile samples, e.g. silicone oil)
- Never operate the hot plate/stirrer in wet areas!
- Be aware of the danger of electric shocks!
-  **Warning!** Be aware of the potential danger of a fire outbreak due to overheating!



-  **Warning!** Wear your personal protective equipment in accordance with the hazard category of the media to be processed. Otherwise there is a risk from:
  - Splashing and evaporation of liquids
  - Ejection of parts
  - Release of toxic or combustible gases.
-  **Warning!** When in an emergency, disconnect the main power plug.
- Gradually increase the speed.
- Reduce the speed if:
  - The medium splashes out of the vessel because the speed is too high
  - The appliance is not running smoothly
  - The container moves on the base plate.
- The safe temperature limit must always be set to at least 25°C lower than the fire point of the media used.
- Beware of hazards due to:
  - Flammable materials
  - Combustible media with a low boiling temperature
  - Glass breakage
  - Incorrect container size
  - Overfilling of media
  - Unsafe condition of container.
- The base plate can heat up due to the action of the drive magnets at high motor speeds, even if the Infrared hot plate / stirrer is not operational.
- Process pathogenic materials only in closed vessels under a suitable extractor hood.
- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- Please observe the operating instructions for any accessories used.
- Ensure that the external temperature sensor (PT 100 ) is inserted in the media to a depth of at least 30 mm.
- The PT 100 external temperature sensor must always be inserted in the media when connected.
- Safe operation is only guaranteed with the accessories described in the "Accessories" chapter.
- Accessories must be securely attached to the device and cannot come off by themselves. The center of gravity of the assembly must lie within the surface on which it is set up.
- When using PTFE-coated magnetic bars, the following has to be noted: Chemical reactions of PTFE occur in contact with molten or solute alkali metals and alkaline earth metals, as well as with fine powders of metals in groups 2 and 3 of the periodic system at temperatures above 300-400°C. Only elementary fluorine, chlorine fluoride and alkali metals attack it; halogenated hydrocarbons have a reversible swelling effect.

### 3.3. For protection of the equipment

- You have received a product designed for industrial and experimental use. Nevertheless, avoid strikes to the housing, vibrations, damage to the operating-element panel, and contamination.
- Make sure that the mains power supply has low impedance to avoid any negative effects on instruments being operated on the same mains.
- Do not expose the unit to sunlight
- Sudden drops may cause damage in the interior of the instrument.
- Transport the instrument with care.
- The device can be damaged when sucking in aggressive gases or vapor through the installed ventilator.
- Press the power button to interrupt the stirrer, rather than disconnect the main power plug directly.
- The voltage stated on the nameplate must correspond to the mains voltage.
- Do not cover the device, even partially e.g. with metallic plates or film. This results in overheating.
- Protect the appliance and accessories from bumps and impacts.
- Ensure that the base plate is kept clean
- Observe the minimum distances between devices, between the device and the wall and above the assembly (min. 800 mm).

## 4. Operating Procedures

### 4.1. Environmental Operating Conditions

The Infrared hot plate / stirrer must operate in the following conditions:

- Indoors
- Altitudes up to 2000 meters
- Temperatures from +5°C to +40°C
- Maximum relative humidity 80% for temperatures up to +31°C, linear decrease down to 50% relative humidity at a temperature of +40°C
- Max. mains fluctuation of  $\pm 10\%$  are permissible
- Protection class according to EN 60 529: IP20
- The unit corresponds to Class I
- Overvoltage category II

## 4.2. Installation

### 4.2.1. Installing the infrared Hot Plate/Stirrer

Place the hot plate/ stirrer on a stable, flat surface and proper environment for operation.

- If a PT100 temperature sensor package was ordered:



Install the sensor holder into the screw inlet on the back of the Infrared hot plate / stirrer



Place the clamp onto the holder and tighten the screw connecting the holder and the clamp




Mount the Pt100 temperature sensor on clamp and tighten the screw connecting the temperature sensor and the clamp



Connect the PT100 temperature sensor cable to the corresponding connector in the back of the Infrared hot plate / stirrer

- 
- Connect the stirrer to the power supply. The power supply voltage, frequency and current are respectively AC 108V~125V、60Hz、10A.
  - Connect the power supply to a power socket with earthing contact.

#### 4.2.2. The dimensions of the infrared hot plates/stirrer (after connect PT100Temperature sensor)




#### Caution!



- Do not use voltages that are higher or lower than 10% of the voltage specified on the label, which is on the backside of the instrument.
- Keep the power cord and temperature sensor cable off of the hot plate while heating.
- Put the solution on the top plate before operating the instrument.
- Heating corrosive liquids under poor ventilation hoods will shorten the life of the electronic components inside the instrument.
- Upon the first heating operation, a particular smell and white smoke can appear. This is normal. Put the instrument under a fume hood and moderately heat for about one hour until the smell and smoke fully disappear.
- If toxic gases are released, air circulation must be kept.
- The safe temperature limit must always be set to at least 25 °C lower than the fire point of the media used.


## 4.3. Operation

### 4.3.1. Overview of the Infrared hot plate / stirrer



No.	Description
1	PT100 Temp. Sensor
2	Holder for PT100Temp. Sensor
3	Heating Top Plate
4	Right Control Knob
5	LCD Digital Display
6	Left Control Knob
7	Instrument Housing
8	Clamp for Temp. Sensor
9	Installation rod for Temp. Sensor

### 4.3.2. Indicators and Functional Elements



No	Icon	Description
1		<b>Left side Bar graph</b> Indicates the control activity of the stirring function,
2		<b>Stirrer symbol</b> Indicates that the stirrer is ready for operation.
3		<b>Right side Bar graph</b> Indicates the control activity of the heating function,
4		<b>Heater symbol</b> Indicates that the stirrer is ready for operation.
5		<b>Residual-heat indicator</b> Warning sign to inform the user that the heating zone is still hot.

### Caution!




- Do not use voltages that are higher or lower than 10% of the voltage specified on the label, which is on the backside of the instrument.
- Keep the power cord and temperature sensor cable off of the hot plate while heating.
- Put the solution on the top plate before operating the instrument.
- Heating corrosive liquids under poor ventilation hoods will shorten the life of the electronic components inside the instrument.
- Upon the first heating operation, a particular smell and white smoke can appear. This is normal. Put the instrument under a fume hood and moderately heat for about one hour until the smell and smoke fully disappear.

### 4.3.3. Operation of the Infrared hot plate / stirrer without a Temperature Sensor


#### 1. Switching the Infrared hot plate / stirrer on



*e.g. the software version*




*Display during self-test*



*e.g. The initial heating stage is always "0"*

- Plug in power
- Switch the Infrared hot plate / stirrer on by pressing the main switch on the right side of the Infrared hot plate stirrer
- The stirrer performs a self-test which is indicated on the display followed by information on the software version
- Press and hold the Right Control Knob for approximately 2 seconds until the Infrared hot plate / stirrer Indicator appears
- The display shows the heating stage "0" (zero).
- Select your desired heating stage within 30 seconds
- If no heating stage is selected, the Infrared hot plate / stirrer switches off again after 30 seconds

## 2. Selecting a Heating Stage




*e.g. the heating stage is "24"*

- If the Infrared hot plate / stirrers operated without a temperature sensor, the heating power will be controlled.
- Select the desired heating stage using the Right Control Knob (turn clockwise to set a higher heating stage turn anti-clockwise to set a lower heating stage). There are 24 heating stages available.
- The Infrared hot plate / stirrer starts immediately after a heating stage other than "0" is selected
- The selected heating stage is shown on the display (e.g. Set to 24 stage), after complete the setting the display will show the Infrared hot plate / stirrer symbols.
- The bar graph indicates the heating activity of the heater.
- The heating element heats up until the selected temperature level is achieved.
- If no heating stage is selected, the Infrared hot plate / stirrer switches off again after 30 seconds
- After having operated in the heating stages 19-24 for 2 hours, the Infrared hot plate / stirrer switches back to the heating stage 18 as a safety function




### 3. Setting the Stirring Speed



*e.g. The initial RPM is always "0"*



*e.g. set RPM to 800*



*e.g. Heating and stirring display simultaneously*

- Press and hold the Left Control Knob for approximately 2 seconds until the Stirrer Indicator appears
- The display shows the initial stirring speed is always at "0" (zero), when switching the unit on.
- The stirring speed can be set from 50 to 1500 rpm. The stirring function can be used alone or together with the heating function.
- By using the left Control Knob to set the stirring speed (turn clockwise to increase the set value, turn anti-clockwise to decrease the set value)
- The bar graph indicates the stirring activity

### 4. Switching the Infrared hot plate / stirrer off



- Press and hold the two Control Knob respectively for approximately 2 seconds until the Heating and stirring Indicator disappears.
- The Infrared hot plate / stirrer and stirrer is now switched off.
- The Residual Heat Indicator continues to light up as long as the glass-ceramics heating zone is still hot
- The installed ventilator continues to operate until the heating zone has cooled down completely.



#### CAUTION!

- Residual heat! Do not touch the heating zone!
- Risk of overheating! Do not pull out the mains plug!
- Do not unplug and turn off the mains of the hot plate before the heating zone has completely cooled down, then turn off the mains switch and pull out the mains plug.

### 4.3.4. Operation of the Infrared hot plate / stirrer with a Temperature Sensor




#### Note!

In contrast to operation without temperature sensor, the laboratory hot plate now features:

- Automatic temperature control instead of fixed heating stages controlled by the heating power.
- Temperature display alternating between set temperature and actual temperature instead of showing the heating stage

#### 1.Connecting the Temperature Sensor



- Make sure that the Infrared hot plate / stirrer is completely switched off
- Be sure to use the correct temperature sensor
- Connect the temperature sensor at the rear of the Infrared hot plate / stirrer
- Make sure that the cable of the temperature sensor is routed so that it cannot touch the heating zone.
- Immerse the temperature sensor into the liquid min.30 mm in depth.


#### 2.Switching the Infrared hot plate / stirrer on



*e.g. The solution temperature is "28"*

- Refer to 4.3.3 for how to set the temperature, except for the following one points.
- Now the instrument will enter the external temperature control mode (to measure and control the temperature of the solution accurately). The maximum set temperature is 300°C.

### 3.Setting the Temperature



- Using the Right Control Knob (turn clockwise to increase the set value, turn anti-clockwise to decrease the set value).



*e.g. The setting temperature is 150°C*



*e.g. The actual temperature is 140°C*

- The bar graph indicates the heating activity of the heater
- The Infrared hot plate / stirrer heats up and maintains the selected temperature
- The display now alternates between the set temperature (with dots between figures) and the actual temperature (no dots between figures) every 5 seconds.

#### 4.Setting the Stirring Speed

Please refer to 4.3.3

#### 5.Switching the Infrared hot plate / stirrer off

Please refer to 4.3.3

### CAUTION!



- This equipment is the source for a strong magnetic field interfering with objects up to a distance of 50cm from the laboratory stirrer. It is therefore recommended to be careful when approaching the stirrer with magneto-sensitive objects such as electronic data carriers (discs, bank cards), mechanical wristwatches or pacemakers, etc.!
- Risk of interference resulting from a magnetic field!
- Keep away magneto-sensitive objects!
- Use a support frame to prevent the stirrer vessel from slipping off the heating zone, if required.

### 4.3.5. Pt100 Temperature Sensor Calibration


The Pt100 temperature sensor can be connected to measure and control the heated liquid temperature. The sensor has been initially calibrated in the factory. If the measured temperature is slightly different from the temperature standard you are using, a follow-up calibration is also possible with the following steps.

#### Calibration Procedure



e.g. the measured temperature "T1" is 24 °C

- Loosen the retaining screw of the instrument and remove the bottom shell (totally 10 pieces) Place the instrument at side view, and remain the switch side.
  - Plug the PT100 temperature sensor into the designated input and turn on the main power switch
  - Put the PT100 sensor and third-party thermometer in the temperature environment of T1, (e.g. 300 ml silicone oil, the temperature stability at 24 degree, this value should be measured and read by the third-party thermometer)
-



e.g. The contact switch on processor board




e.g. display shows "S1" and it's corresponding AD read value



e.g. set "S1" to "T1" ( 24°C)



e.g. display shows "S2" and it's corresponding AD read value



e.g. set "S2" to "T2" ( 151°C)



- Find the processor board in instrument and press the contact switch 5 times to enter the calibration procedure;

- Then LCD screen light up and the right side digital display shows "S1", "S1" means Set 1(The first setting temperature for calibration).
- the left side display shows its corresponding AD read value T1

- Wait for the temperature to reach and equilibrium 24°C. (measured by the third-party thermometer)
- refer to the shows value on third-party thermometer and twist the right side control knob to set the actual value of T1 (e.g. Set to 24 )

- Put the PT100 temperature sensor and third-party thermometer in temperature environment of T2, ( e.g. 300 ml silicone oil, the temperature stability at 151°C , ( this value should be measured and read by the third-party thermometer )
- Press the contact switch again, the display shows "S2" "S2" means Set 2(The second setting temperature for calibration).
- the left side display shows its corresponding AD read value of T2

- Wait for the temperature to reach an equilibrium
- refer to the shows value on third-party thermometer and twist the right side control knob to set the actual value of T2 (e.g. Set to 151°C )

- Press the contact switch to complete the calibration.

## 4.4. The USB interface

The instrument can be connected to a PC through the COM port which is using the USB connection, it can be plugged directly to the PC's available USB connection.

**1. We provide USB port COM port driver and installation instructions, the driver can be downloaded from the following address.**

Driver download address:

<http://www.ftdichip.com/Drivers/VCP.htm>

Installation instructions download address: <http://www.ftdichip.com/Support/Documents/InstallGuides.htm>

### 2. RS-232 command set

Item	Input / output	Content of "#"	Meaning	Remarks
Set the equipment parameters				
Set temperature unit	set to # ↵	0	°C	
Set temperature	out_sp_00 ### ↵	Decimal number	—	Unit: °C or °F (set by equipment)
Set speed	out_sp_01 ### ↵	Decimal number	—	Unit: RPM
Set state	out_mode_05 # ↵	0	stop	
		1	start	
Timing sending once every 3 seconds, Temperature unit: °C, It has nothing to do with the equipment settings				
internal control	rs ### ### ↵	Decimal number	Internal temperature, speed	eg: rs 1000 120 means current Internal temperature 100.0°C, speed 120RPM
External control	prs ### ### ### ↵	Decimal number	External temperature, Internal temperature, speed	eg: prs 800 1000 120 means current External temperature 80.0°C, Internal temperature 100.0°C, speed 120RPM
Send "status" will return the following information				
Current setting temperature	in_sp_00 ### ↵	Decimal number	—	Unit: °C or °F (set by equipment)
Current setting speed	in_sp_01 ### ↵	Decimal number	—	
Current set temperature unit	set to # ↵	0	°C	
		1	°F	
Current setting	status # ↵	0	stop	


status		1	start	
When the alarm will be sent to the following information				
alarm	status # ↙	03	Internal temperature exceeds the upper limit	
		04	Control temperature over-temperature	
		06	Motor speed anomaly	
		07	Safety temperature overrun	

## 5. Cleaning and Maintenance

### 5.1. Routine Cleaning

The device is maintenance-free.

#### Cleaning



For cleaning disconnect the main plug.

Only use cleansing agents which have been recommended by *WIGGENS*

Use to remove:

Dyes	isopropyl alcohol
Construction materials	isopropyl alcohol/water containing surfactant
Cosmetics	isopropyl alcohol/water containing surfactant
Foodstuffs	water containing surfactant
Fuels	water containing surfactant

- Do not allow moisture to get into the appliance when cleaning.

- Wear protective gloves when cleaning the devices.

- Before using another than the recommended method for cleaning or decontamination, the user must ascertain with *WIGGENS* that this method does not destroy the instrument.



#### **Note:**

Do not use chlorine bleach, chlorine-based cleanser, abrasives, ammonia, steel wool or scouring pads with metal content or similar harsh solvents or abrasives. These may damage the surface of the instrument.



## 5.2. Maintenance

Do not attempt to service or repair a *WIGGENS* Infrared hot plate / stirrer. If the Infrared hot plate / stirrer housing is opened the warranty becomes void. Contact *WIGGENS* for return authorization and return instructions.

Ordering spare parts

When ordering spare parts, please give:

- Machine type
- Manufacturing number, see type plate
- Item number and designation of the spare part.

Repair

Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards. For this, use the "certificate of compliance" form which you can obtain from *WIGGENS*. If your appliance requires repair, return it in its original packaging. Storage packaging is not sufficient when sending the device - also use appropriate transport packaging.

## 6. Transport and Storage

- Clean the Infrared hot plate / stirrer so that it is free from any materials which may be harmful to the health. Provide a material safety data sheet where appropriate.
- Place the infrared hot plate / stirrer unit and its parts into the original packing or a container with necessary protection to prevent damage during transport. Seal the original packing or container with packing tape.
- Store the packed unit in a dry place.



### CAUTION!

Failure to clean, maintenance, and handle the infrared hot plate / stirrer as outlined can lead to damages or be harmful to the health.











## 7. Accessories and Spare Parts







### 7.1 Temperature Sensor and Holder

Model	Description	Order No.
PT100 Temperature sensor, Type I	Length: 170 mm; Diameter: 4 mm; Material: Stainless steel; Admissible temperature: -30 ~ +300°C	PT100-01
PT100 Temperature sensor, Type II	Length: 300 mm; Material: Stainless steel	PT100-02
PT100 Temperature sensor, Type III	Length: 150 mm; Material: Stainless steel, PTFE coated	PT100-03
PT100 Temperature sensor, Type IV	Length: 300 mm; Material: Stainless steel, PTFE coated	PT100-04
PT100 Temperature Sensor, Type VI	Length: 250 mm; Diameter: 4 mm; Material: Glass; Admissible temperature: -30 ~ +300°C	PT100-06
Holder for Temperature Sensor	Holder and clamp for PT100 temperature sensor; Suitable for SLR / WH240-PLUS/WH240-HT / WH220-PLUS/WH220-HT	PT100-05

## 7.2. Stir Bar and Retriever Examples

Model	Description	Order No.
Cylindrical Stirrer Bars		001.110.6
Plain Stir Bars		001.210.6
Octahedral Stir Bars		001.513-R/B/Y
Micro Stir Bars		001.802
Colored Micro Stir Bars		001.802-R/B/Y
Oval Stir Bars		001.610
Octaoval Stir Bars		001.3319
Cylindrical Stir Bars With a Removable Pivot Ring		001.1712
Tapered Stir Bars		001.1910
Double Ended Stir Bars (Natural)		001.1335

Model	Description	Order No.
Test Tube Wings		001.2201.1
Cross-Shaped Stir Bars		001.2401
Crosshead Stirrer Magnets(Double Sided)		001.1110
Crosshead Stirrer Magnets(Single Sided)		001.1110.1
Hub Stirrer Magnets		001.2301
Tube Stirrer Magnets		001.1609
Disc-Shaped Stirrer Magnets		001.709
Spherical Stirrer Magnets		001.1512
Glass Encased Stir Bars		001.1206
Extra Large Stir Bars		001.0120

Double Ended Stir Bars (Colored)		001.1335-R/B/Y	Stir Bar Sets (18 Bars)		001.3003
Triangular Stir Bars		001.412	Rare Earth 'Turbo' Stir Bars		001.2610.RE
Ribbed Triangular Stir Bars		001.1812	Stir Bar Retrievers		004.150



For more information about Accessories please contact your local supplier



### CAUTION!

For safety and guarantee reasons only original accessory parts are to be used!

## 8. Service

### 8.1. Trouble-Shooting

Cause	Remedy
After switching on the unit, the display shows no light and the Infrared hot plate / stirrer does not react to any input.	<ol style="list-style-type: none"> <li>1. Ensure that the mains electricity plug is plugged into a working socket outlet and check if the main switch is in the "on" position.</li> <li>2. Open the fuse holder of the power cord, which you can find at the back of the instrument. If the fuse is damaged, replace it with a 4A / 230V fuse. Clean the holder before your replacement.</li> <li>3. If the fuse is not damaged and the malfunction cannot be determined, please contact the <i>WIGGENS</i> support.</li> </ol>
After switching on the unit, the power switch is lit up, but the display shows a blank screen.	<ol style="list-style-type: none"> <li>1. This is probably a malfunction of the control board. Please contact the <i>WIGGENS</i> support.</li> </ol>



*WIGGENS* reserves the right to carry out technical modifications with repairs for providing improved performance of the instrument.

### 8.2. Warranty

In accordance with *WIGGENS* warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs. The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

*WIGGENS* reserves the right to decide the validity of any warranty claim. In case of faults arising either due to faulty materials or workmanship, parts will be repaired or replaced free of charge.

Any other compensation claims, such as consumables, damages caused by corrosion or accidental breakage, are excluded from this guarantee.

This warranty may only be altered by a specifically published amendment. No individual has authorization to alter the provisions of this warranty policy or its amendments.

## 8.3. Contact /Technical Service

If your device is not working properly:

⇒ Please inform *WIGGENS* Instruments by using our contact information.

You have contacted *WIGGENS* Instruments?

⇒ Copy and complete the Confirmation of condition of unit from these operating instructions.

⇒ Please repack the device appropriately for transport and send to *WIGGENS* Instruments together with the Confirmation of condition of unit.

### Our contact details

#### **WIGGENS GmbH**

Add: Gässlesweg 22-24, 75334 Straubenhardt, Germany

Tel.: 0049 7248 4529088

#### **Wiggins Co., Ltd.**

Room 426, Hall A, Office Building M8, No.1 Jiuxianqiao East Road, Chaoyang District, Beijing 100015, China

Tel: +86 400-809-2068

Fax: +86 400-809-2068-112

info@wiggins.com

service@wiggins.com

www.wiggins.com

## Confirmation of condition of unit

In the case of repair, copy and complete the Confirmation of condition of unit and send it to *WIGGENS* Instruments.

### 1. Details about the unit

Product number

---

Serial number

---

Reason for repair

---



---

### 2. Has the device been cleaned, decontaminated/sterilized?

Yes

No

### 3. Is the unit in a condition which does not represent any health threats for the staff of our service department?

Yes

No

If not, which substances has the unit come into contact with?

---

### 4. Legally binding declaration

The customer is aware of being legally liable to *WIGGENS* Instruments for any damages arising from incomplete and incorrect information.

Date

Signature

Company stamp

---

## Please Note

The shipper is responsible for the return of the goods in well-packed condition, suitable for the mode of transport.

### Sender information

Name

---

Company

---

Department, research group Street

---

Zip code, city

---

Country

---

Phone

---

E-mail

---



**WIGGENS GmbH**

Gässlesweg 22-24, 75334 Straubenhardt, Germany

Tel.: 0049 7248 4529088

**WIGGENS China**

Room 426, Hall A, Office Building M8, No.1 Jiuxianqiao East Road,

Chaoyang District, Beijing 100015, China

Tel: +86 400-809-2068

Fax: +86 400-809-2068-112

[info@wiggens.com](mailto:info@wiggens.com)

[www.wiggens.com](http://www.wiggens.com)