

TM Series Vortex Shakers TM-2F



- Very accurate shaking speed due to micro-chip speed controller
- Stable and reliable performance
- High shaking speed of up to 2500 rpm
- Shaking operation can be switched between continuous and on-touch operation
- Rotation direction can be switched
- Suitable for shaking test tubes or other small volume vials
- Large range of platforms and other accessories
- Aluminum housing
- Suitable for deep well titer plates
- Can be used as a stirrer and flask shaker

Specifications

Model	TM-2F
Shaking Speed(rpm)	0-2500
Max. Volume (L) in Stirring Mode	2
Housing Material	Plastic-Aluminum
Shaking Amplitude (mm)	4.5
Speed Control/Display	Turning Knob LED Display
Operation Mode	Continuous / On-Touch
Accessories	Accessory Pack
Power (W)	17.3
Dimension (mm)	115X230X60
Order No.	36012F

Accessories for TM Series

Universal Tray

Can be used with different foam plastic platforms for simultaneous shaking of several test tubes

Order No.: 3601-01



Universal Base Plate

Base plate for foam plastic platforms with 120 mm diameter

Order No.: 3602-01



Magnetic Stirring Plate

Can be directly installed on the shaker and used for stirring of smaller samples

Order No.: 3601-02



Test Tube Platform (54 Holes)

Can be used in combination with the universal base plate for shaking operation of 54 Eppendorf tubes

Order No.: 3602-02



Elisa Plate Tray

Can be used with an ELISA Plate

Order No.: 3601-03



Test Tube Platform (18 Holes)

Can be used in combination with the universal base plate for shaking operation of 18 test tubes with 10 mm diameter

Order No.: 3602-03



Test Tube Platform (6 Holes)

Can be used in combination with the universal tray for shaking operation of 6 test tubes with 12 mm diameter

Order No.: 3601-04



Test Tube Platform (12 Holes)

Can be used in combination with the universal base plate for shaking operation of 12 test tubes with 12 mm diameter

Order No.: 3602-04



Test Tube Platform (4 Holes)

Can be used in combination with the universal tray for shaking operation of 4 test tubes with 16 mm diameter

Order No.: 3601-05



Test Tube Platform (8 Holes)

Can be used in combination with the universal base plate for shaking operation of 8 test tubes with 16 mm diameter

Order No.: 3602-05



One-hand insert

80mm

Order No.: 3601-06



Test Tube Platform (8 Holes)

Can be used in combination with the universal base plate for shaking operation of 8 test tubes with 20 mm diameter

Order No.: 3602-06



Flask Holder

Can be used in combination with the universal base plate for shaking operation of Erlenmeyer or round flasks with a volume between 100 and 250 ml

Order No.: 3602-07



WIGGENS GmbH

Gässlesweg 22-24, 75334 Straubenhardt, Germany

Tel.: 0049 7248 4529088

WIGGENS China

Room 303, Hall C, Office Building M8, No.1 Jiuxianqiao East
Road, Chaoyang District, Beijing 100015, China

Tel: +86 400-809-2068

Fax: +86 400-809-2068-112

Email: info@wiggens.com service@wiggens.com

Website: www.wiggens.com